
3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Eosin Y disodium salt - Dye content ≥85 %

Product Number:	**E6003**
CAS Number:	17372-87-1
MDL:	MFCD00005040
Formula:	C20H6Br4Na2O5
Formula Weight:	691.85 g/mol



TEST	Specification
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Solubility (Color)	Orange to Very Dark Orange and Red to Very Dark Red and Orange-Red and Red-Orange
1 mg/ml H2O	
Solubility (Turbidity)	Clear
Dye Content	≥ 85 %

Specification: PRD.1.ZQ5.10000034863

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.