Product Specification Sheet

VWR DOMED PCR CAP STRIPS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.


VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 732-1518

Material: Polypropylene

Dimensions :



ISO 9001
ISO 14001
OHSAS 18001

Product Specification Sheet

VWR DOMED PCR CAP STRIPS

Format: Strips of 8 caps

Compatible with: 96-well PCR plates as well as standard and low profile strip tubes

Colour: Transparent

Quality: Free from DNase, RNase og DNA

Storage: In closed original cartons, dry, at room temperature

Weight: Approx. 1.0 g / piece

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **October 2012**

ISO 9001
ISO 14001
OHSAS 18001