Technical Datasheet

VWR® CLEANROOM POLYESTER WIPER 2PLY SB 23X23CM

EU cat. No 114-9323

Intended use

- Cleaning the critical working surface in clean room

Characteristics

- Constructed from continuous-filament polyester in interlock knitting.
- Ultrasonic Bonded 2-layer construction.
- Processed and packed in class 10 cleanroom.
- With Border-sealed exhibits lowest particle and fibers levels than other cutting technologies.
- Higher absorbency compared to 2 single wipers.

General Specifications

Item		Test method	Standard
Material		-	100% Polyester (75 Denier)
Weight (g/m^2)		Balance	220 ± 10
Dimension (cm)		Ruler	23 x 23
Liquid Particle Count (>=0.5um,counts/m^2)		IEST-RP-CC004.3	$\leqslant 10.0x10^7$
Non Volatile Residue-NVR (g/m^2)	DI-Water	IEST-RP-CC004.3	$\leqslant$ 0.02
	IPA	IEST-RP-CC004.3	$\leqslant$ 0.04
Absorbency (ml/m^2)		IEST-RP-CC003.4	$\geqslant$ 700
Ions Contamination (ppm)	Cl-	IEST-RP-CC004.3	$\leqslant$ 0.50
	Na+	IEST-RP-CC004.3	$\leqslant$ 0.50
	K+	IEST-RP-CC004.3	$\leqslant$ 0.50

Packaging and Storage

Package/Quantity	100pcs/poly bag, 10bags/case
Best storage Conditions	Room temperature, dry and shady
Shelf life	5 years
Country of Origin	China