





Bottles Dimension

Catlog #	Description	Capacity (ml)	Bottle Shape	Bottle Material	Closure Material	Approx. Brim Volume (ml)	Colour	Closure Size, mm	Torque (in-lb)		Height with Closure	Height without closure	Size (L) (Length X Width)	Neck ID	Product Weight	Min Wall Thickness	Leakproof	Autoclavable
									Min	Max	mm	mm	mm	mm	Gms	mm		
215-4550	WIDE MOUTH SQUARE BOTTLE HDPE 60ml	60ml	Square	HDPE	PP	70ml	Natural	28-415	16	21.5	83	80.9	35.5 X 35.5	20.8	13.5	0.7	Yes	-
215-4551	WIDE MOUTH SQUARE BOTTLE HDPE 175ml	175ml	Square	HDPE	PP	185ml	Natural	38-415	22	30	103.5	100.5	51 X 51	28.5	32	1	Yes	-
215-4552	WIDE MOUTH SQUARE BOTTLE HDPE 250ml	250ml	Square	HDPE	PP	295	Natural	43-415	28	38	114.5	111.5	60 X 60	32.7	35	0.8	Yes	-
215-4553	WIDE MOUTH SQUARE BOTTLE HDPE 500ml	500ml	Square	HDPE	PP	575	Natural	53-415	33.5	44.5	144	141.5	72.7 X 72.7	43.7	61.5	1	Yes	-
215-4554	WIDE MOUTH SQUARE BOTTLE HDPE 1000ml	1000ml	Square	HDPE	PP	1150	Natural	63-415	50.2	55.2	181	177	90.5 X 90.5	52.9	106	1.1	Yes	-
215-4555	WIDE MOUTH SQUARE BOTTLE PP 60ml	60ml	Square	PP	PP	70ml	Natural	28-415	16	21.5	83	81.5	35.5 X 35.5	21	13.5	0.7	Yes	Yes
215-4556	WIDE MOUTH SQUARE BOTTLE PP 175ml	175ml	Square	PP	PP	185ml	Natural	38-415	22	30	103.5	100.5	52 X 52	28.7	32	0.8	Yes	Yes
215-4557	WIDE MOUTH SQUARE BOTTLE PP 250ml	250ml	Square	PP	PP	295	Natural	43-415	28	38	114.5	112.5	60 X 60	33	35	0.8	Yes	Yes
215-4558	WIDE MOUTH SQUARE BOTTLE PP 500ml	500ml	Square	PP	PP	575	Natural	53-415	33.5	44.5	144	141.5	75 X 75	43.7	61.5	1	Yes	Yes
215-4559	WIDE MOUTH SQUARE BOTTLE PP 1000ml	1000ml	Square	PP	PP	1150	Natural	63-415	50.2	55.2	181	178	93.5 X 93.5	53.2	105.7	1.2	Yes	Yes
215-4560	WIDE MOUTH SQUARE BOTTLE HDPE AMBER 60	60ml	Square	HDPE	PP	70ml	AMBER	28-415	16	21.5	83	80.9	35.5 X 35.5	20.8	13.5	0.7	Yes	-
215-4561	WIDE MOUTH SQUARE BOTTLE HDPE AMBER 17	175ml	Square	HDPE	PP	185ml	AMBER	38-415	22	30	103.5	100.5	51 X 51	28.5	32	1	Yes	-
215-4562	WIDE MOUTH SQUARE BOTTLE HDPE AMBER 25	250ml	Square	HDPE	PP	295	AMBER	43-415	28	38	114.5	111.5	60 X 60	32.7	35	0.8	Yes	-
215-4563	WIDE MOUTH SQUARE BOTTLE HDPE AMBER 500ml	500ml	Square	HDPE	PP	575	AMBER	53-415	33.5	44.5	144	141.5	72.7 X 72.7	43.7	61.5	1	Yes	-
215-4564	WIDE MOUTH SQUARE BOTTLE HDPE AMBER 1000ml	1000ml	Square	HDPE	PP	1150	AMBER	63-415	50.2	55.2	181	177	90.2 X 90.2	52.9	107.4	1.1	Yes	-

Material

Bottle : Refer above table for material type, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

Closure : Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Bottle Application Temperature:

- HDPE: -100° C to +120° C
- PP : -40° C to +135° C

Autoclave Condition

- Temperature : 121°C
- Time : 20 min
- Pressure: 15 psi

Autoclave Caution:

1. Before Autoclaving separate / disengage the caps threads from bottle.
2. Do not stack products, this may cause vacuum and may lead to physical damage / collapse of product.

Material Recycle Code



Description: Wide Neck Square Bottle

Drawing: Not to Scale

avantor™

Date: 02/15/2023

Revision: 0