

ThermoFisher SCIENTIFIC

Product Specifications

Product **397691000**
Ethanol,99.5%,Extra Dry,absolute, AcroSeal®

General Product Data

Version	0
CAS No	64-17-5
Molecular weight	46.07
Molecular formula	C2 H6 O
Linear formula	CH3CH2OH
Flash point (°C)	12

Product Specifications

Appearance (Color)	Clear colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms
Water	=<0.005 %
Residue after evaporation	=<0.001 %
GC	>=99.5 %

Product Specifications are subject to amendment and may change over time.

Acros Organics BV

ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.thermofisher.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329