SP SCIENCEWARE
Bel-Art · H-B Instrument · Wilmad-LabGlass
661 Route 23 South
Wayne, New Jersey 07470
Tel: 973-694-0500 * Fax: 973-694-7199
Internet Address: http://www.bel-art.com

Technical Data Sheet

Product #: 14650-0000

Description: Urbanti High-Speed Filter Funnel

Material: TPX = Polymethylpentene (TPX RT 18XB)

Dimensions: Top ID 196mm, Stem ID 20mm, Stem OD 26mm, Stem Length 160mm.

Other: Steam Autoclavable at 121°C for 15 minutes
Materials meet FDA 21CFR177.1520, USP Class VI
Manufactured: Bel-Art Products Pocomoke City, MD USA