Technical sheet

Saffron in Alcoholic Solution

Code product: **369200**

Specifications:

• **Appearance:**	Clear orangey solution
• **Absorption Maximum in 95°-Ethanol:**	433 nm ≤ λ ≤ 434 nm

Applications:

Histo-Cytology: used as a specific staining solution for collagen in Haemalum-Eosin-Saffron and Haemalum-Phloxin-Saffron techniques for routine trichromic stainings of histological sections. Alcoholic Saffron stains collagen golden yellow.

Mycology: used in P.A.S. staining for the detection of some protozoa and the parietal structures of fungi, *mycelium* and yeasts in smears (cutaneous squama, oropharynx and œsophagus cells, broncho-alveolar lavage liquid).

Updated: 25/04/2013

RAL Diagnostics - Site Montesquieu - 33651 MARTILLAC Cedex – France - Tel : 33+557 960 404 - Fax : 33+557 960 405 - www.ral-diagnostics.fr