Product Specification / Certificate of Conformity

Product Code	BLH0150P
Product Description	Azlon® Wide neck bottle, HDPE
Closure Material	Polypropylene
Material particular	HPDE
Closure Colour	Natural
Colour	Translucent
Capacity ml	150
Pack quantity	10

Dimensions and Tolerances

Diameter mm	51
Neck ID mm	32
Height mm	116

Product Compliance

Latex	Yes
Reach	Yes
TSE_BSE Compliant	Yes
Food Compliant	EU 10/2011 / FDA CFR 21 177

This is to certify that the goods detailed above have been inspected and tested in accordance with our BS EN ISO9001 and BS EN ISO13485 quality systems, and conform in all respects to the acceptable quality levels (AQL) detailed therein and any relevant specifications or drawings, unless recorded below as concessionary goods.

Kevin Marshall
Head of Quality

Date Issued: 20/09/2021

DWK Life Sciences Limited
Unit 4, Riverside, 2 Campbell Road,
Stoke-on-Trent, ST4 4RJ,
UK

Phone: +44 (0)1782 444406
Fax: +44 (0)1782 940436
certificates@dwk.com
www.dwkltd.com