Type in Product Names, Product Numbers, or CAS Numbers to see suggestions.

Specification Sheet

Product Name	Ricinoleic acid analytical standard
Product Number	90969
Product Brand	SIAL
CAS Number	141-22-0
Molecular Weight	298.46

TEST	SPECIFICATION
APPEARANCE (COLOR)	Colorless to Faint Yellow
APPEARANCE (FORM)	Liquid
PURITY (GC AREA %)	≥ 99.0 %
SPECIFIC ROTATION (20/D)	5.0 - 8.0 Degrees
CONCENTRATION	C=2 in ACETONE
REFRACTIVE INDEX N20/D	1.469 - 1.473
PROTON NMR SPECTRUM	CONFORMS TO STRUCTURE
EXPIRY PERIOD	36 MONTHS

EXPIRY PERIOD	36 MONTHS
TEST	SPECIFICATION

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
Reproduction of any materials from the site is strictly forbidden without permission.
Site Use Terms | Privacy Policy | General Terms and Conditions of Sale | Copyright Consent | Cookies Settings