Seit 1927
über 90 Jahre

Since 1927
more than 90 Years

Amarell
GmbH & Co. KG
Präzisions-Thermometer / Aräometer

Lindenstraße 43
D-97892 Kreuzwertheim

Tel. +49 (0) 9342 9283-0
Fax: +49 (0) 9342 39860

Internet: www.amarell.de
E-Mail: thermometer@amarell.de

Technical Data Sheet

Art. No.:	A300730-CB/WC
Type:	Thermometer similar ASTM 54C, solid stem, white backed
Specification:	Congealing Point
Measurement range:	20+100,6°C
Subdivision:	0,2°C
Length:	310mm
Diameter:	6-7mm
Filling:	Capillary tube specially coated inside, with non-wetting blue special liquid
Immersion:	Total
Certificate:	Works Certificate with 4 test points

September'2021