Technical Data Sheet

Article number	**44302**
Description	**Wooden applicators**

Art.-No.	Length	Applicator Ø	Cotton End
44302010	approx. 15 cm	approx. 2.2 mm	approx. 12 x 5 mm Ø
44302020	approx. 20 cm	approx. 2.2 mm	approx. 15 x 5 mm Ø
44302030	approx. 23 cm	approx. 2.2 mm	approx. 15 x 5 mm Ø

Properties: with cotton end

Tel.: +49 9779/808-0
Fax: +49 9779/808-88
Stettener Str. 22-24
97647 Sondheim/Rhön · Germany
info@assistent.eu
www.assistent.eu

Glaswarenfabrik Karl Hecht GmbH & Co KG | Amtsgericht Schweinfurt HRA 8063 | USt.IdNr./VAT-No.: DE132198535
Persönlich haftende Gesellschafterin: Hecht GmbH | AG Schweinfurt HRB1160
Geschäftsführer: Karl-Heinz Hecht, Thomas Diepholtz | Einzelprokura: Tonja Müller