
InLab® Expert Pro

pH Ref °C ISM

● ● ● ○

For general applications in the lab and outdoor;
robust and maintenance-free



Sensor type	combined pH electrode
Order number	51343101

Specifications

pH range	0...14
zero point (pH)	7.00 ± 0.25
temperature range	0...100 °C
temperature probe	NTC 30 kΩ
shaft material	PEEK
membrane glass	U
membrane resistance	< 250 MΩ
reference system	ARGENTHAL™
reference diaphragm	double open junction
reference electrolyte	XEROLYT® polymer
bridge electrolyte	-
electrode head	-
cable length	1.2 m
connector	BNC / RCA (Cinch)
cell constant (cm^{-1})	-
cell type	-
indicating metal	-

Storage

InLab storage solution (30111142)
or FRISCOLYT-B® (51340053)

© 12/2015 MTII www.mt.com/electrode-guide

METTLER TOLEDO