FREUDENBERG
MEDICAL


Technical Drawing Tube Standard Silicone (TSS)

Wacker Elastosil 4305/60

Materialnumber	OD (B)	±	ID (A)	±	Wall (C)	Length
00 45 63 41 40- (60-411-40)	0,64 mm	0,05 mm	0,31 mm	0,05 mm	0,17 mm	15 m
00 45 63 41 41- (60-411-41)	0,94 mm	0,05 mm	0,51 mm	0,05 mm	0,22 mm	15 m
00 45 63 41 42- (60-411-42)	1,19 mm	0,05 mm	0,64 mm	0,05 mm	0,28 mm	15 m
00 45 63 41 43- (60-411-43)	1,65 mm	0,05 mm	0,76 mm	0,05 mm	0,45 mm	15 m
00 45 63 41 44- (60-411-44)	2,16 mm	0,08 mm	1,02 mm	0,08 mm	0,57 mm	15 m
00 49 30 05 29- (60-411-45)	1,96 mm	0,08 mm	1,47 mm	0,08 mm	0,25 mm	15 m
00 45 63 41 46- (60-411-46)	2,41 mm	0,08 mm	1,58 mm	0,08 mm	0,42 mm	15 m
00 45 63 41 47- (60-411-47)	3,18 mm	0,08 mm	1,58 mm	0,08 mm	0,80 mm	15 m
00 49 30 05 30- (60-411-48)	3,18 mm	0,08 mm	1,98 mm	0,08 mm	0,60 mm	15 m
00 45 63 41 49- (60-411-49)	4,88 mm	0,10 mm	2,64 mm	0,10 mm	1,12 mm	15 m
00 45 63 41 50- (60-411-50)	4,65 mm	0,10 mm	3,35 mm	0,10 mm	0,65 mm	15 m
00 45 63 41 51- (60-411-51)	6,35 mm	0,10 mm	3,18 mm	0,10 mm	1,59 mm	15 m
00 45 63 41 52- (60-411-52)	7,95 mm	0,18 mm	4,78 mm	0,18 mm	1,59 mm	15 m
00 45 63 41 53- (60-411-53)	9,53 mm	0,18 mm	4,78 mm	0,18 mm	2,38 mm	15 m
00 45 63 41 54- (60-411-54)	9,53 mm	0,18 mm	6,35 mm	0,18 mm	1,59 mm	15 m
00 45 63 41 55- (60-411-55)	11,13 mm	0,25 mm	6,35 mm	0,18 mm	2,39 mm	15 m
00 45 63 41 56- (60-411-56)	12,70 mm	0,25 mm	6,35 mm	0,18 mm	3,18 mm	15 m
00 45 63 41 57- (60-411-57)	12,70 mm	0,38 mm	7,94 mm	0,25 mm	2,38 mm	15 m
00 45 63 41 58- (60-411-58)	12,70 mm	0,38 mm	9,53 mm	0,25 mm	1,59 mm	15 m
00 45 63 41 59- (60-411-59)	14,30 mm	0,38 mm	9,53 mm	0,25 mm	2,39 mm	15 m
00 45 63 41 60- (60-411-60)	15,88 mm	0,38 mm	9,53 mm	0,25 mm	3,18 mm	15 m
00 49 30 05 31- (60-411-61)	17,48 mm	0,38 mm	12,70 mm	0,25 mm	2,39 mm	15 m
00 45 63 41 62- (60-411-62)	19,05 mm	0,38 mm	12,70 mm	0,25 mm	3,18 mm	15 m
00 45 63 41 63- (60-411-63)	22,23 mm	0,38 mm	15,88 mm	0,25 mm	3,18 mm	15 m
00 49 30 05 32- (60-411-64)	23,83 mm	0,38 mm	15,88 mm	0,25 mm	3,98 mm	15 m

Wacker Elastosil 4305/80

Materialnumber	OD (B)	±	ID (A)	±	Wall (C)	Length
00 49 30 20 73- (60-831-01)	9,00 mm	0,25 mm	6,00 mm	0,25 mm	1,50 mm	15 m
00 49 30 18 72- (60-831-02)	14,00 mm	0,51 mm	10,00 mm	0,38 mm	2,00 mm	15 m
00 49 30 12 70- (60-831-03)	25,00 mm	0,51 mm	19,00 mm	0,51 mm	3,00 mm	15 m
00 49 30 09 45- (60-831-05)	16,65 mm	0,51 mm	9,53 mm	0,38 mm	3,56 mm	15 m



Freudenberg Medical Europe GmbH
Liebigstr. 2-8
67661 Kaiserslautern
Telefon +49 (0)631 5341-7500
Fax +49 (0)631 5341-7600
info@freudenbergmedical.de

10.10.2017,

Quality Assurance Freudenberg Medical Europe GmbH, Sven Thiel