Product Specification Sheet

VWR BOTTLES, WIDE MOUTH, WITH SCREW CAP, PP, TRANSLUCENT

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity (ml)	Ø ext.(mm)	Height (mm)	Neck I-Ø (mm)	Neck Ø ext. (mm)	Pk
215-5685	1000	96	206	50	65	12
215-5682	100	47	94	24	32	10
215-5683	250	63	119	31	40	10
215-5684	500	78	152	39	50	12
215-5681	50	39	77	24	32	10

Material : Polypropylene, Translucent

Density(at 20°C): 0.89 - 0.91g/cm³

Other Characteristics:

- Perfect for sample collection, storage and transport. Ideal for liquids, pastes and granular material.
- The wide opening permits easy filling.
- Can be filled with hot substances. With easy to read blue graduations.
- Very good chemical resistance
- Heat resistant up to +121 °C (for 20 minutes), prolonged use +100 °C
- Leakproof

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt, ask us.*

Date of issue: **August 2014**