
A part of: ThermoFisher SCIENTIFIC

National Scientific
197 Cardiff Valley Rd
Rockwood, TN 37854

Technical Data Sheet

Product Number: 6AK40AOTAS
Description: EPA Vial and Closure Kit

Component	***Specification***
Vial Physical Properties	
Material	51 expansion clear borosilicate glass
Height	95mm
Diameter	27.5mm
Thread Finish	24-400
Wall Weight	1.1mm
Cap/Septa Physical Properties	
Cap material	polypropylene
Cap size	24-414
Type	open top screw cap
Cap Color	white
Septa Material	PTFE/silicone
Rubber	silicone
Rubber Color	clear
Film	PTFE
Film Color	natural
Septa thickness	2.54mm
Septa diameter	22.7mm

	9/15/23	
Revision Number	0	
Approved By	Patrick Giltnane	