HIRSCHMANN®

Technical data sheet

Article: Ceramus classic 1 – 5 ml

Code no.: 933 00 00

Standard: ISO 8655

Tolerances:

Accuracy A:	± 0,6%
Coefficient of variation CV:	0,2%

Material:

Piston:	highly pure Al_2O_3 – ceramic (99,7%)
Cylinder:	Borosilicat - glass 3.3
Valve block:	Halar (ECTFE)
Valve ball:	Borosilicate glass 3.3
Valve seat:	Borosilicate glass 3.3
Valve spring:	Hastelloy

Eberstadt, 04.05.2021

.....
Roland Benz/ Quality Manager