Motic®

MORE THAN MICROSCOPY



ST-30C 2LOO

This Classic Motic Series is specially designed and developed for both the learning environment as well as quick quality control management. You can choose between the innovative cordless LED illumination available over the standard tungsten system making your images brighter and sharper without producing heat.

Model	ST-30C-2LOO
Optical System	Greenough
Observation Tube	Binocular head
Interpupillary distance	51-79mm
Diopter adjustment	On the left eyepiece, +/- 5 diopter
Inclination	45° inclined
Eyepieces	WF10X/20mm
Objectives system	Turret changer
Objectives magnification	2X-4X
Working distance	47mm
Stand	Pole type stand with transmitted illumination
Base	190x115mm
Column/Arm	130mm high
Head holder	For Ø13mm column with incident illumination
Focus mechanism	Coarse focusing system with tension adjustment
Focusing stoke	42mm
Incident illumination	Tungsten 12V/10W
Transmitted illumination	Tungsten 12V/10W
Transformer	Internal
Power supply	220V
Accessories included	Black, white and frosted stage plates, dust cover
Length x Width x Height	190x115x285mm
Weight	2,7Kg

GO DIGITAL

Add a Moticam camera to capture, document, annotate and share images and videos with the Motic Images Plus software.

Please verify the local availability of the product with your Motic Representative.