S MURRAY & CO
Holborn House Old Woking
Surrey GU22 9LB England

Telephone +44 (0) 1483 740099
Fax +44 (0) 1483 755111/724462
e-mail sales@smurray.co.uk

www.smurray.co.uk

SPECIFICATION SHEET – CO60UC20

Bottle Extra Wide Mouth 60ml With R3/51mm Cap - Assembled

Product Code CO60UC20

Dimensions (Bottle only)

Volume/Capacity	60ml
Height	52.4mm (+/- 0.8mm)
Diameter	51.6mm (+/- 0.8mm)

Neck Finish

Neck type	Screw neck, R3/51mm
Cap reference	G808/11
Bottle reference	G806/03

Material

Bottle	Soda Glass, Clear
Cap	Urea, PVDC Wad 51mm/R3

Packing Specification

Pack quantity	20 units
Pack weight	1.9kg
Pack dimensions	28.5 x 21 x 6cms

Product Image



52.4mm
(+/- 0.8mm)
(not incl. cap)

51.6mm (+/- 0.8mm)

Product Image – Complete Pack of 20



This specification sheet is the property of S Murray & Co Ltd

Signed: Geoff Hambly Dip RSA, Quality Director

Dated 1st January 2018