Technical Data Sheet

Article number	43600
Description	Fridge-freezer thermometer

43600010	
Measuring range	-40 °C up to +50 °C
Dimensions	30 x 20 x 210 mm
Weight	25 g
Material	Plastic
Steril	No
Packaging	1 Piece

Advantages: - rotatable

- with cooling zones
- High accuracy due to 2-point adjustment

Execution: The cylindrical thermometer is equipped with a hook for hanging.

Glaswarenfabrik Karl Hecht
GmbH & Co KG
Stettener Str. 22-24
97647 Sondheim/Rhoen

Tel: +49 9779 8080
Fax: +49 9779 808 88
info@assistent.eu
www.assistent.eu

USt. ID: USt.IdNr./VAT-No.:
DE 132198535
Amtsgericht Schweinfurt HRA 8063
Persönlich haftende Gesellschafterin:
Hecht GmbH AG Schweinfurt HRB1160

Geschäftsführer:
Steffen Wachter
Thomas Diepholtz
Einzelprokura:
Tonja Müller