Technical Data Sheet

Article number	**42779070**
Description	**Flat bottom tubes**
Material	Soda glass
Capacity	18.5 ml
OD x Height	22 x 70 mm
Wall thickness	1 mm
Bottom	Flat bottom
Rim	Smooth rim
Pack of (Pcs)	100

Our conformity certified flat bottom test tubes are commonly used in laboratories and are designed to stand upright without the need for additional support. Their stable base facilitates easy handling and prevents tipping over during experiments or sample storage. We offer these tubes in two different sizes, catering to diverse scientific needs.

Glaswarenfabrik Karl Hecht
GmbH & Co KG
Stettener Str. 22-24
97647 Sondheim/Rhoen

Tel: +49 9779 8080
Fax: +49 9779 808 88
info@assistent.eu
www.assistent.eu

USt. ID: USt.IdNr./VAT-No.:
DE 132198535
Amtsgericht Schweinfurt HRA 8063
Persönlich haftende Gesellschafterin

Geschäftsführer:
Steffen Wachter
Thomas Diepholtz
Prokura: Tonja Müller