


REV	DESCRIPTION OF CHANGE	DRAWN BY	DATE	ECN
C	REVISE "E" DIMENSION	BL	01/29/10	20100012009513228
–	REDRAWN ON DWK TEMPLATE AND STANDARDS	EL	07/29/20	N/A







FINISH DETAIL
SCALE 2:1

NOTES:
1. BOTTLE CAPACITY: 68.0 ±2.5cc OVERFLOW.
2. BOTTLE WEIGHT: 9.0 ±0.8g.

ITEM	QTY.	MATERIAL: DESCRIPTION AND SPECIFICATIONS
2	72	NATURAL POLYPROPYLENE
1	72	NATURAL POLYPROPYLENE

THIRD ANGLE PROJECTION

TITLE: 60 ML STAR BTL, NM, PP NAT, 20–410 CAP

THIS CONFIDENTIAL AND PROPRIETARY DRAWING AND THE INFORMATION IN IT ARE THE EXCLUSIVE PROPERTY OF DWK LIFE SCIENCES. IT IS PROVIDED WITH THE UNDERSTANDING AND INTENTION THAT NOTHING CONTAINED HEREIN WILL BE REPRODUCED OR DISCLOSED TO THIRD PARTIES WITHOUT THE PRIOR EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES.

ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

TOLERANCES
DIMENSIONS IN INCHES UNLESS SPECIFIED
.XX ± .01"
.XXX ± .005"
.XXXX ± .0005"
FRACTIONAL ± 1/64"

DURAN WHEATON KIMBLE

ITEM NUMBER: 209166	DRAWING NUMBER: CSD–209166	DRAWN BY: BL	DRAWN DATE: 04/04/08		
SAP NUMBER: 209166	OLD MATERIAL NUMBER: N/A	SHEET SIZE: B	SHEET NUMBER: 1 OF 1	SHEET SCALE: 1:1	