Certificate of Conformity

VWR Catalogue Number	612-2703
Description	VWR®, Safety Bulb Pipette Fillers
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-08-04

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared that the products are inspected and controlled through whole-production processing in accordance with current applicable product specifications and quality control SOPs.

Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. inspected in according to AQL Level II 4.0 sampling plan per MIL-STD-1916.

Product Specifications

Material	Bulb is injection molded of vulcanized rubber; hard adapter is injection molded of acrylonitrile-butadiene-styrene (ABS) with non-toxic white pigment; soft adapter is injection molded or silicone. Elastic tie.
Sterilization	Non-Sterile.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	These products are not made from natural rubber latex.

BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	According to EU Directive 2015/863 No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Temperatures: 15 °C to 38°C (59 °F to 100 °F) Relative humidity during storage (non-condensing: 10% to 90%)
Shelf-Life	Indefinite shelf life.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	08-2025