CHAIR LAB STOP AND GO CASTORS 560-810MM



ECN	**132-1377**	
Dimensions in mm	Backrest:	H 420
	Seat:	H 560-810 W450 D 420
Packing data	Weight:	12,0 kg
	Volume:	0,19 m³

Basic equipment		**Benefit**
Laboratory compliance	Washable und resistant to disinfectant. GMP and cleanroom compatible. Can be used in laboratories with safety classification S1 to S3. Minimal joints and hermitec (Hygienic-Design).	Suitable for laboratories and all hygiene-critical environments particularly health care, universities, biotechnology, chemistry and pharmaceuticals.
Handling	Practical handle on backrest. Extremely lightweight. Compact size.	Ensures maximal mobility and ease of use.
Backrest	Synthetic backrest with 3D-Flex function.	The backrest adapts to the movements of the user and therefore follows the upper body when leaning back or turning sideways. This ensures a dynamic seating position and provides optimal support for the body in any position. Intuitive and automatic.
Seat	Flex function at front edge of seat.	The front edge of the seat automatically follows the forward-angled working position of the body. This ensures optimal blood circulation and effectively prevents pressure pain.
Seat Shell	Plastic seat shell in blue	Attractive design. Improves the ambience of the workplace and employee satisfaction.
Upholstered Seat	Magic artificial leather: High-quality synthetic leather, washable,low-maintenance and impervious to oil and disinfecting agents. Colour blue.	The soft and comfortable upholstery offers a very high sitting comfort. However, damages caused by sharp parts can not be excluded. Ideal for the use in production and laboratory..
Seat Height	Seat height infinitely adjustable between 560-810 mm	Large variability of seat height adjustment.
Base	Made of polyamide black with foot ring	Extremely high stability, durability and static stability.
Castors	Soft Stop & Go-castors for hard floors	Offer a superb level of stability. Blocks the chair when seated. Easily movable when there is no pressure on the chair.
Foot ring	Height adjustable, chrome-plated foot ring.	Wide foot plate provides good support. Height can be quickly and easily adjusted.
Standards	Suitable for cleanrooms category 3 in accordance with DIN EN ISO 14644-1 GMP-Conform Laboratory safety classification S1 to S3 Tested for harmful substances technical standard DIN 68 877-1 GS-Certified for highest levels of safety Blue Angel environmental certificate Sanitized Quality management DIN EN ISO 9001 Environmental Management DIN EN ISO 14001 Verified Environmental Management EMAS Occupational Health & Safety OHSAS 18001 EcoVadis CSR Rating 2016 Gold.	Verified by certificates from independent testing institutes, quality, environmental friendliness and procurement security.

Guarantee	10-year long-term-guarantee	Free shipping for spare parts during the whole period of validity

Labsit