Product Technical Data & Conformity Report

Date: 18/03/2025
Product Type: Bottom Weld Bag
Product Ref: w/o 277658 o/n 75102
Application: Hazardous Waste

1. Manufactured Size

Type	Closed Width (MM)	Open Width (MM)	Length (MM)	Thickness MU
B/W Bag		790	1020	125

2. Dimension, Properties & Others

Item	Standard	Method
Closed Width	0	Tape Measure
Open Width	790	Tape Measure
Length	1020	Tape Measure
Thickness MU	125	Dial Gauge &/Or Weight Per Metre
Blend	VMED S 30	Blend Sheet
Blend Properties	LDPE Based Med/High Slip With 30%+ Regran	Blend Sheet
Masterbatch	Yellow1300977	Batch Labels
Recycling Category	LDPE 4	Batch Labels
Plastic Tax	Exempt From Plastic Tax	Legislation
Print Design	JENCONS BIOHAZ X 25 (121-029)	Laysheet
Ink Colour(s)	1 colour black	Pantone Book
Quantity Per Box	75	Machine Index
Packaging	Bags Boxed Pallet(s) Top Sheeted & Stretch Wrapped	Production &Warehouse Operations