Laboratory glassware,
apparatus and
precision instruments

Glaswarenfabrik Karl Hecht GmbH & Co KG

Stettener Str. 22 - 24
97647 Sondheim v. d. Rhön
Germany

Tel. +49 9779 808-0
Fax +49 9779 808-88

Email: info@assistent.eu
Internet: http://www.assistent.eu
VAT No. DE 132 198 535

Data sheet / Quality certificate

Date of issue: 19.12.2019

Product description

Brushes
galvanized steel wire, with loop or wooden handle

natural bristles: boiled, undyed - water absorption about 40% - resistant up to 150 °C
nylon bristles: resistant up to 100 °C

Country of origin: Germany HS code 96039099

Product code	ex	**Total length** mm	**Brush** length mm	diam. mm	**Bristles**	**Tip**	
91120000000	8000	900	150	13	natural	protective cap	for burette
91120000001	8001	150	50	8	natural	cotton end	for test tube
91120000002	8002	220	50	15	natural	cotton end	for test tube
91120000003	8003	300	90	30	natural	cotton end	for test tube
91120000004	8004	400	80	4	natural		for pipette
91120000029	8004A	400	90	3	natural		for pipette
91120000005	8005	250	70	8	natural		for pipette
91120000006	8006	400	100	12	natural		for pipette
91120000007	8007	250	70	15	natural	bent	for pipette
91120000008	8008	500	110	8	natural		for pipette
91120000009	8009	850	110	8	natural		for pipette
91120000010	8010	300	90	30	natural	bent	
91120000011	8011	300	110	50	natural	bent	
91120000012	8012	260	60	25	natural		for butyrometer
91120000013	8013	430	130	80	natural / nylon		for beaker
91120000014	8014	390	125	85/60	natural / nylon		wooden handle for beaker

Laboratory glassware,
apparatus and
precision instruments

**Glaswarenfabrik
Karl Hecht
GmbH & Co KG**

Stettener Str. 22 - 24
97647 Sondheim v. d. Rhön
Germany

Tel. +49 9779 808-0
Fax +49 9779 808-88


Email: info@hecht-assistent.de
Internet: http://www.assistent.eu
VAT No. DE 132 198 535

Data sheet / Quality certificate

Date of issue: 19.12.2019

Product code	ex	**Total length** mm	**Brush** length mm	diam. mm	**Bristles**	**Tip**	
91120000017	8017	450	140	40	natural	tied tip	for bottle
91120000018	8018	500	150	55	natural	tied tip	for bottle
91120000019	8019	600	160	80	natural	tied tip	for bottle
91120000020	8020	500	140	100	natural	tied tip	for bottle
91120000021	8021	300	65	11	nylon	bent	
91120000022	8022	300	65	16	nylon	bent	
91120000023	8023	300	65	20	nylon	bent	
91120000024	8024	300	70	28	nylon	bent	
91120000025	8025	300	75	38	nylon	bent	
91120000026	8026	300	100	75	nylon	bent	
91120000027	8027	1000	100	25	natural	tied tip	for burette
91120000028	8028	800	120	18	natural	tied tip	for burette

Photos (91120000002) (91120000018)





Glaswarenfabrik
Karl Hecht
GmbH & Co KG