
Technical Data Sheet



VWR® Exact Volume Transfer Pipettes

Description: Exact Volume Pipet


Purpose: Transfer liquid

Materials: LDPE DJ200A

Features

- Uses: under -20 ~ 50 °c
- Sterilization: No
- BSE/TSE Statement: Meet the requirements of BSE/TSE.
- Products are manufactured under the ISO 13485 standard.

VWR® Transfer pipettes



EU Part No.	Graduations (ml)	Scale type	Total Length (mm) ± 2mm	Exact Volume （ul） ± 20ul	Bulb draw (ml)	Drop/ml ±1	individually Pack	Packing
612-5846	N/A	N/A	83	100	N/A	30	NO	500/PK, 5000/Case
612-5847	N/A	N/A	102	250	N/A	22	NO	500/PK, 5000/Case
612-5848	N/A	N/A	115	200	N/A	23	NO	500/PK, 5000/Case
612-5849	N/A	N/A	59	50	N/A	28	NO	500/PK, 5000/Case
612-5850	N/A	N/A	75	75	N/A	28	NO	500/PK, 5000/Case