Sodium Hydroxide, 1.0 Normal
Volumetric Solution

Material No.: 4693-09
Revision No.: 0

Product Specification

Test	Specification
Normality	0.995 - 1.005

For Laboratory,Research,or Manufacturing Use
Protect from air to avoid absorption of carbon dioxide.

Packaging Site: Paris Mfg Ctr & DC

Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Mansford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700