Ethyl Acetate
ChromAR HPLC

Material No: H078-25
Effective from: 2013/07/02

Product Specification

Test	Specification
Assay (GC)	>= 99.8 %
Water (KF)	<= 0.05 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.0005 %
Ultraviolet Absorbance (1.00-cm cell vs. water):	
260 nm	<= 0.15
270 nm	<= 0.05
280 nm	<= 0.02
300 nm	<= 0.004
UV Cut-off	<= 252 nm

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000