Material No.: 5616-02
Revision No.: 0

Product Specification

Test	Specification
Appearance (clear, colorless solution)	Passes Test
Color (APHA)	≤ 10
Normality	1.995 – 2.005
Sulfate (SO_4)	≤ 1 ppm
Ammonium (NH_4)	≤ 3 ppm
Trace Impurities – Heavy Metals	≤ 0.5 ppm
Trace Impurities – Iron (Fe)	≤ 0.1 ppm

For Laboratory,Research,or Manufacturing Use
Standardization at 25°C traceable to NIST Standard Reference Material.

Packaging Site: Paris Mfg Ctr & DC

James Ethier
Jamie Ethier
Vice President Global Quality