**Reference:** 499263OO **Technical Data Sheet**

**Product:** **STAINLESS CONNECTOR - MEDIA IN BAGS**

## Specification

Connector suitable for culture media in bags.

## Presentation

| | Packaging Details | Shelf Life | Storage |
|---|---|---|---|
| 1 unit | 1 Connector + inner box + louter box +label. | 122 months | 1-40 ºC |

with: ± 5 ml

## Composition

Stainless steel connectors for ready-bags.

## Description /Technique

The culture bag has a connection point with a penetrable cap (injection port) latex-free polycarbonate, for any additive injection required, and for connecting the bag to any standard equipment laboratory dosing with the connector. The injection of the connector must be performed aseptically after its previous sterilization, and then connected to a synthetic rubber silicone tube with 4.8 mm(0.63 inch) of internal diameter.

## Quality control

### Physical/Chemical control

Color : Greyish pH: at 25ºC

### Microbiological control

Not Performed - Only Sterility Test

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020

Not Applicable

| Microorganism | Growth |
|---|---|
| *pseu* | NOT PROCEED |

### Sterility Control

Not Performed - Chemical Reagent without nutritive properties.

## Bibliography

N. A.