
Technical Product Data

Technical Data Sheet | Labcon Product Number 3651-875-000-9



Dimensions are milimeters. Do not scale



Image is not actual size. Colors may differ.

SuperClear® 1.5 mL Freestanding Screw Cap Microcentrifuge Tubes with Caps, Sterile

Fluid Measure	mL	**Freezable**	Yes
Fluid Volume	1.5	**Freeze Temp**	-80°C
Height (mm)	45.4	**Autoclavable**	121°C 15min
Width (mm)	10.4	**Sterile**	Yes
Length (mm)		**Sterile Processing**	Certified as sterile with a sterility assurance level (SAL) of 10-6
Pack Qty	500	**Centrifugation**	Yes
Case Qty	5000	**Centrifugation Speed (RCF)**	20000
Package Type	Resealable Bags	**Shelf Life (Months)**	36
Material	Polypropylene	**Certification**	
Cap Material	Polypropylene	**Regulatory**	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is **© Copyright 23-Mar-2022** . All Rights Reserved. Technical Data Sheet Version 1/2020