Ethylene Glycol
Baker Analyzed

Material No: 7037.9025
Effective from: 2013/06/10

Product Specification

Test	Specification
Assay (by GC)	>= 99.0 %
Acidity (as CH3COOH)	<= 0.01 %
Boiling Point	195-198°C
Boiling Range	194.0-200.0°C
Color (APHA)	<= 10
Density (g/ml) at 20°C	1.110-1.114
Residue after Ignition	<= 0.005 %
Water (H2O)	<= 0.2 %
Trace Impurities (in ppm)	
Chloride (Cl)	<= 5
Iron (Fe)	<= 0.2

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000