G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

A Geno Technology, Inc. (USA) brand name

Adenosylhomocysteine Nucleosidase


A Recombinant Enzyme

(Cat. # 786-427)



INTRODUCTION

Adenosylhomocysteine nucleosidase (EC 3.2.2.9) converts S-adenosylhomocysteine (AdoHcy) to S-(deoxy-D-ribos-5-yl)-L-homocysteine + adenine. A highly purified recombinant enzyme protein supplied at a 10µM concentration




ITEM(S) SUPPLIED (Cat. # 786-427)

Description	Size
Adenosylhomocysteine nucleosidase (EC 3.2.2.9) (10µM)	0.5ml

RELATED PRODUCTS

Download our Bioassays Handbook.



http://info.gbiosciences.com/complete-bioassay-handbook

For other related products, visit our website at www.GBiosciences.com or contact us.

Last saved: 8/10/2012 CMH



www.GBiosciences.com