







Dichtung Ø79± 0,3 x 3,1±0,3
Art.-Nr. 11 529 20 (Lose)
Art.-Nr. 11 529 21 (5er-P)

X





Trade mark Ø2,8

Schrift: Arial. 5mm hoch, 20pt
0,3mm erhaben

Schrift: Arial. 4mm hoch, 16pt
0,3mm erhaben

Entformungsschräge 0,75°

Weitergabe sowie Vervielfältigung dieses Dokuments, Verwertung und Mitteilung seines Inhalts sind verboten, soweit nicht ausdrücklich gestattet. Zuwiderhandlungen verpflichten zu Schadenersatz. Alle Rechte für den Fall der Patent-, Gebrauchsmuster- oder Geschmacksmustereintragung vorbehalten.
The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

Index	Aenderung / changes	Datum / date	Name
E	Art.-Nr. 1165889 has been removed	26-FEB-2020	Schneider
D	Change of logo to DURAN	28-NOV-2017	Schneider
C	Art.-Nr 11 658 89 ohne Dichtung hinzu	16-SEP-2014	Weindorf
B	Artikel der Form angepasst	22-MAI-09	Weindorf

zul. Abweichung allowed dev.	ISO 2768-m K		
Oberfläche surface	ISO 1302		
	Datum / date	Name	
bearb. drawn	22.05.2009	Weindorf	
geprüft checked	29.11.2017	A. Rees	

DURAN WHEATON KIMBLE
DURAN Produktions GmbH & Co. KG
www.DWK-LifeSciences.com

Werkstoff / material	Stück quantity	Gewicht weight	Pos. position
PSU-Compound		~55g	

Benennung / title: Schnellschraubverschlusskappe GLS80 — Maßstab scale: 4:1

Art.-Nr. 11 658 88, 11 664 16 (mit Dichtung)

Form-Nr. mould number: GLS80 — A3

Zeichnungsnummer / drawing number: 0743104-GZA-GLS80 — Revision: E — Blatt/sheet: 1 /1

Ers. / repl.f: 0743104-GZA-GLS80-C — Ers. / repl.d:

Schutzvermerk ISO 16016 beachten / Refer to protection notice ISO 16016 — INV