Acetonitrile

ChromAR HPLC

Material No: 6711-25
Effective from: 2013/03/06

Product Specification

Test	Specification
Assay	>= 99.8 %
Water (KF)	<= 0.03 %
Titrable Acid (meq/g)	<= 0.0005
Titrable Base (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.0005 %
Ultraviolet Absorbance (1.00-cm cell vs. water):	
190 nm	<= 0.52
200 nm	<= 0.05
210 nm	<= 0.03
220 nm	<= 0.02
230 nm	<= 0.01
240 nm	<= 0.005

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000