IKA

designed for scientists



SY stand

/// Data Sheet

The EasySyn reactor systems with double-walled reactor vessels from 500 ml to 5000 ml provide high flexibility using a single, specially designed stand system.

- Quick and easy exchange of different reactor vessels
- Simple alignment of the overhead stirrer
- Adjustment of the operating height to the individual body size.

IKA

designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	500 x 950 x 580
Weight [kg]	20
Permissible ambient temperature [°C]	5 - 4 0