Product Specification Sheet

VWR SPREADER, T-SHAPED

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release.The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Colour	Packed	Pk
612-2651	T-shaped spreader	Blue	1/peel pouch	500
612-2652	T-shaped spreader	Blue	5/peel pouch	1.000
612-2653	T-shaped spreader	Blue	10/peel pouch	1.000

Material: Plastic, Sterile

Other characteristics:

- Designed for spreading and dispersing liquids across the surface of agar culture plates.
- Used for performing bacterial (for example, CFU) plate counts on water, milk and other liquid samples.
- Smooth rounded surfaces to prevent cutting or gouging of agar
- T-shaped to ensure even pressure is applied while spreading

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **April 2013**