# Heidelberger Extension Tubings



**Description of used materials (simple):**

Heidelberger Extension Tubing (4097300)

| | **Designation** | **Polymer** | **Polymeric Material** |
|---|---|---|---|
| 1. | Luer lock connector with protective cap | | |
| 1.1. | Protective cap | PE-LD | Polyethylene, low density |
| 1.1. | Luer lock fitting male | MABS | Methyl methacrylate-acrylonitrile-butadiene-styrene plastic |
| 2. | Luer lock connector | | |
| 2.1. | Luer lock fitting male | MABS | Methyl methacrylate-acrylonitrile-butadiene-styrene plastic |
| 3. | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| 3. (alternative) | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| | | PVC-P | Poly(vinyl chloride) plasticized |
| 4. | Luer lock fitting female | PVC | Poly(vinyl chloride) |
| 5. | Luer lock adapter | PVC | Poly(vinyl chloride) |
| 6. | Closing cone Luer-Lock | PE-HD | Polyethylene, high density |
| 7. | Closing cone | PE-HD | Polyethylene, high density |

**Description of used materials (simple):**

Heidelberger Extension Tubing (4097173)

| | **Designation** | **Polymer** | **Polymeric Material** |
|---|---|---|---|
| 1. | Luer lock connector with protective cap | | |
| 1.1. | Protective cap | PE-LD | Polyethylene, low density |
| 1.2. | Luer lock fitting male | MABS | Methyl methacrylate-acrylonitrile-butadiene-styrene plastic |
| 2. | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| 2. (alternative) | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| | | PVC-P | Poly(vinyl chloride) plasticized |
| 3. | Luer lock adapter | PVC | Poly(vinyl chloride) |
| 4. | Closing cone | PE-HD | Polyethylene, high density |

# Heidelberger Extension Tubings



**Description of used materials (simple):**

Heidelberger Extension Tubing (4097408)

| | **Designation** | **Polymer** | **Polymeric Material** |
|---|---|---|---|
| 1. | Luer lock connector with protective cap | | |
| 1.1. | Protective cap | PE-LD | Polyethylene, low density |
| 1.2. | Luer lock fitting male | MABS | Methyl methacrylate-acrylonitrile-butadiene-styrene plastic |
| 2. | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| 2. (alternative) | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| | | PVC-P | Poly(vinyl chloride) plasticized |
| 3. | Luer lock adapter | PVC | Poly(vinyl chloride) |
| 4. | Closing cone | PE-HD | Polyethylene, high density |

# Heidelberger Extension Tubings



**Description of used materials (simple):**

Heidelberger Extension Tubing (4097262)

| | **Designation** | **Polymer** | **Polymeric Material** |
|---|---|---|---|
| 1. | Luer lock connector with protective cap | | |
| 1.1. | Protective cap | PE-LD | Polyethylene, low density |
| 1.2. | Luer lock fitting male | MABS | Methyl methacrylate-acrylonitrile-butadiene-styrene plastic |
| 2. | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| 3. | Luer lock adapter | PVC | Poly(vinyl chloride) |
| 4. | Closing cone | PE-HD | Polyethylene, high density |

# Heidelberger Extension Tubings



**Description of used materials (simple):**

Heidelberger Extension Tubing (4097190)

| | **Designation** | **Polymer** | **Polymeric Material** |
|---|---|---|---|
| 1 | Connector female with LL | | |
| 1.1 | Luer lock fitting female | PVC | Poly(vinyl chloride) |
| 1.2 | Closing cone Luer-Lock | PE-HD | Polyethylene, high density |
| 2 | Tubing | PVC-P | Poly(vinyl chloride) plasticized |
| 3 | Luer lock connector with protective cap | | |
| 3.1 | Protective cap | PE-LD | Polyethylene, low density |
| 3.2 | Luer lock fitting male | MABS | Methyl methacrylate-acrylonitrile-butadiene-styrene plastic |