KORFF AG
Niedermattstrasse 35
CH-4538 Oberbipp
Telefon +41 (0)32 636 33 32
Fax +41 (0)32 636 23 09
info@korff.ch
www.korff.ch

Product specification for art. no. 60278

Product description: Aluminium foils for laboratories, soft, plain

Product specification:

Material:	pure aluminium foil, soft,
Thickness:	0,030 mm
Sheet dimension:	110 x 120 mm, with interleaf paper
Colour/Design:	silver / silver
Art. No.:	60278



(only reference picture, not art. specific)

Packing: 1000 pcs/unit
approx. 1,4 kg per packet (indicated weight is including interleaf paper)

Product qualification: Aluminium is stable between the pH values 4.5 and 8.5.
Do not use this product for storage or processing of acidic, basic or salty foods. However, solubility increases markedly at pH below 4.5 and above 8.5 (Elinder and Sjögren, 1986).

Heat resistant: Melting point at 580 °C - 600 °C (for sheets not coloured / silver)

Storage and durability: Temperature 5 – 30 °C, max. 70% relative humidity, no direct sun exposure

All data are based on our knowledge and experience. It's intended as an advice without legally binding.

Oberbipp, 12. January 2023