Technical sheet

Harris Haematoxylin

Code product: **361070**

Specifications:

• **Appearance:**	Purplish blue red solution
• **Absorption Maximum in water:**	530 nm ≤ λ ≤ 560 nm

Applications:

Hormonal Cytology:

- used as a nuclear staining solution in Papanicolaou staining, regressive variation for the trichromic staining of vaginal smears. Harris Haematoxylin colours nuclei violet-blue more or less dark.
- used as a nuclear staining solution in Harris-Shorr staining, regressive variation for the bichromic staining of vaginal smears. Harris Haematoxylin colours nuclei violet-blue or black.

Cytology : used in simplified Harris-Shorr staining for spermocytograms (morphological analysis of spermatozoa).

Note:

Filter the solution before use.

Updated: 25/04/2013

RAL Diagnostics - Site Montesquieu - 33651 MARTILLAC Cedex – France - Tel : 33+557 960 404 - Fax : 33+557 960 405 - www.ral-diagnostics.fr