
IKA

designed for scientists



AS 260.2 Fixing clip attachment

/// Data Sheet

For shaking flasks, Erlenmeyer flasks and bottles with a round cross-section (without fixing clips).
Dimensions (W x D x H): 330 x 330 x 24 mm
Capacity: 56 x AS 2.1 (25 ml), 23 x AS 2.2 (50 ml), 23 x AS 2.3 (100 ml), 11 x AS 2.4 (250 ml), 9 x AS 2.5 (500 ml), 5 x AS 2.6 (1000 ml)
Weight: 1.290 g

Clamping screws for fixing the attachment are included in the scope of delivery of the shaker (HS / KS 260, KS 3000).

IKA

designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	330 x 24 x 330
Weight [kg]	1.294
Permissible ambient temperature [°C]	5 - 40