

Original HDK Catheter Clamp

HDK no.	**741000, 741100, 743000, 741200**

HDK no.	Description	Material
741000	Standard Packaging	POM – polyoxymethylene The product is latex free
741200	Minigrip	
743000	Bulk Packaging	

Description	Catheter Clamp for catheters/tubes. **NB:** Using the Catheter Clamp as an umbilical cord clamp is prohibited.	
Dimensions	58 x 13 mm (closed)	
Quality control	Visual inspection Mil. std. 105E Inspection Level S4 Shrinkage AQL 0.65 Closing test AQL 1.00 Burrs AQL 4.00 Impurities AQL 4.00	
Sterilization	The clamps can be autoclaved (steam sterilized) or sterilized with EtO. Warning: The clamps do not withstand irradiation.	
CE marked	Medical Device, Class I	
Storage	Must be dry at temperatures (-10° – +30°C)	
Packaging	Standard 200 pcs. per box – 2,400 pcs. per export carton Minigrip 100 pcs. per bag – 3,000 pcs. per export carton Bulk 3000 pcs. in a polybag per export carton The packaging materials are latex free	
Manufacturer	Knudsen Plast A/S, Denmark	ISO 9001 and ISO 13485 certified
for	H. Dam Kærgaard A/S, Denmark	ISO 13485 and ISO 14001 certified
Version date	1. marts 2017	
	Subject to changes	