Drop-Boy



Bottle with drop dispenser and separate screw cap.

- Cut off dropper tip before first use
- Reclosable
- Precise drop dosing

Bottle carrier



For 6 bottles up to max. 95 mm Ø.

- PP
- Suitable for food

Bottle carrier box



For 8 bottles up to max. 95 mm Ø. With closed floor which can be used as spill trough.

- HDPE
- Inside measurements: Length 330 mm, width 100 mm
- Collecting volume: 3.9 litres per side

Drop-Boy

Capacity ml	Dia. mm	Height capless mm	Item no.	
10	24	50	0307-0010	N3
25	31.3	64.5	0307-0025	N3
50	40	77	0307-0050	N3

Bottle carrier

No. of bottles	Container dia. mm	Outer dimension of rim L x W x H mm	Item no.	
6	95	310x210x310	0361-0006	N3

Bottle carrier box

No. of bottles	Container dia. mm	Outer dimension of rim L x W x H mm	Item no.	
8	95	400x290x215	0361-0008	N3