				Eco Side access rack							
Freezer brand	Co-brand	Model/size	Box config (height x depth)	ECN	Rack dimensions (HxDxW in mm)	Number of 2" boxes per rack	Number of compartments/freezer	Number of racks/compartment	Number of racks/freezer	Number of 2" boxes	
VWR Collection	Second generation	471-1228 (300 box)	5X5	479-2276	280x691x140	25	4	3	12	300	
VWR Collection	Second generation	471-1229 (400 box)	5X5	479-2276	280x691x140	25	4	4	16	400	
VWR Collection	Second generation	471-1230 (500box)	5X5	479-2276	280x691x140	25	4	5	20	500	
VWR Collection	Second generation	471-1231 (600 box)	5X5	479-2276	280x691x140	25	4	6	24	600	
VWR Collection	Eco Premium	Eco Premium 352	5x4	479-2275	280x562x140	20	2	4	8	160	352
VWR Collection	Eco Premium		6x4	479-2277	321x562x140	24	2	4	8	192	
VWR Collection	Eco Premium	Eco Premium 528	5x4	479-2275	280x562x140	20	2	6	12	240	528
VWR Collection	Eco Premium		6x4	479-2277	321x562x140	24	2	6	12	288	
Thermo	HERA	HDE30086V	4x5	479-2274	225x691x140	20	5	3	15	300	
Thermo	HERA	HDE40086V	4x5	479-2274	225x691x140	20	5	4	20	400	
Thermo	HERA	HDE50086V	4x5	479-2274	225x691x140	20	5	5	25	500	
Thermo	HERA	HDE60086V	4x5	479-2274	225x691x140	20	5	6	30	600	
Thermo	REVCO	RDE30086V	4x5	479-2274	225x691x140	20	5	3	15	300	
Thermo	REVCO	RDE40086V	4x5	479-2274	225x691x140	20	5	4	20	400	
Thermo	REVCO	RDE50086V	4x5	479-2274	225x691x140	20	5	5	25	500	
Thermo	REVCO	RDE60086V	4x5	479-2274	225x691x140	20	5	6	30	600	
Thermo	FORMA	FDE30086V	5X5	479-2276	280x691x140	25	4	3	12	300	
Thermo	FORMA	FDE40086V	5X5	479-2276	280x691x140	25	4	4	16	400	
Thermo	FORMA	FDE50086V	5X5	479-2276	280x691x140	25	4	5	20	500	
Thermo	FORMA	FDE60086V	5X5	479-2276	280x691x140	25	4	6	24	600	
Thermo	Thermo Scientific	TDE30086V	5X5	479-2276	280x691x140	25	4	3	12	300	
Thermo	Thermo Scientific	TDE40086V	5X5	479-2276	280x691x140	25	4	4	16	400	
Thermo	Thermo Scientific	TDE50086V	5X5	479-2276	280x691x140	25	4	5	20	500	
Thermo	Thermo Scientific	TDE60086V	5X5	479-2276	280x691x140	25	4	6	24	600	
Thermo	TSX	TSX400	5X5	479-2276	280x691x140	25	4	4	16	400	
Thermo	TSX	TSX500	5X5	479-2276	280x691x140	25	4	5	20	500	
Thermo	TSX	TSX600	5X5	479-2276	280x691x140	25	4	6	24	600	
Thermo	TSX	TSX700	5X5	479-2276	280x691x140	25	4	7	28	700	
Eppendorf (New Brunswick)	CryoCube	F570h	4x4	479-2273	225x562x140	16	5	5	25	400	
Eppendorf (New Brunswick)	CryoCube	F740i 3 comp	7x4	479-2278	391x562x140	28	3	6	18	504	
Eppendorf (New Brunswick)	CryoCube	F740i 5 comp	4x4	479-2273	225x562x140	16	4	6	24	384	528
Eppendorf (New Brunswick)	CryoCube		6x4	479-2277	321x562x140	24	1	6	6	144	
PHCBi (Panasonic)	TwinGuard	MDF-DU502VX	6x4	479-2277	321x562x140	24	4	4	16	384	
PHCBi (Panasonic)	TwinGuard	MDF-DU702VX	6x4	479-2277	321x562x140	24	4	6	24	576	
Stirling Ultracool	TwinGuard	SU780XLE	4x5	479-2274	225x691x140	20	6	5	30	600	
Stirling Ultracool	TwinGuard	SU105EU	4x3	479-2272	225x427x140	12	2	3	6	72	
Binder	UF	UF V 500	5x4	479-2275	280x562x140	20	2	4	8	160	352
Binder	UF		6x4	479-2277	321x562x140	24	2	4	8	192	
Binder	UF	UF V 700	5x4	479-2275	280x562x140	20	2	6	12	240	528
Binder	UF		6x4	479-2277	321x562x140	24	2	6	12	288	