Swann-Morton® LIMITED

Penn Works
Owlerton Green
Sheffield, S6 2BJ
Telephone: (0114) 234 4231
Email: info@swann-morton.com
Web: www.swann-morton.com

STERILE DISPOSABLE SCALPEL WITH GUARD SPECIFICATIONS

Supplied in boxes of 10

Blades in this range:

10, 10A, 10R, 11, 12, 15, 20, 21, 22, 23, 24

COMPONENT	MATERIAL
BLADE	Stainless Steel **ANALYSIS:** CARBON: 0.60% - 0.70% SILICON: 0.50% max MANGANESE: 1% max CHROMIUM: 12% - 13.5% SULPHUR: 0.025% max PHOSPHORUS: 0.03% max NICKEL: 0.5% max
HANDLE	Medium Impact Polystyrene
GUARD	Low Density Polyethylene
PAPER	Medical kraft paper grid lacquered
FILM	Polyamide / Tie Layer / Low Density Polyethylene
UNIT CONTAINER BOX	White lined chipboard with litho machine varnish

All Components are manufactured in Europe.

All packets and boxes are marked with the current lot number and expiry dates for traceability; boxes are bar coded where applicable.

Products are sterilised by Cobalt 60 Gamma Irradiation by:-Swann-Morton (Services) Limited, Owlerton Green, Sheffield S6 2BJ