

ThermoFisher SCIENTIFIC

Product Specifications

Product **201930250**
Eosin Y,pure,Indicator grade

General Product Data

Version	2
CAS No	17372-87-1
Molecular weight	691.85
Molecular formula	C20 H6 Br4 Na2 O5
Linear formula	

Product Specifications

Appearance (Color)	Dark red-brown
Appearance (Form)	Crystalline powder
Infrared spectrum	Conforms
UV spectrum	Complies
Loss on drying	=<10 %
Specific extinction E (1%/1cm)	>=1100 (lambda max, 0.001M Na2CO3) on dried substance
UV Lambda max	511 to 519 nm

Product Specifications are subject to amendment and may change over time.

Acros Organics BV

ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.thermofisher.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329