DWK LIFE SCIENCES
DURAN WHEATON KIMBLE

DWK Life Sciences Limited

Unit 4, Riverside 2, Campbell Road Stoke-on-Trent, ST4 4RJ, UK

Registered in England and Wales No.7723574

Tel: +44 (0)1782 444406
Fax: +44 (0)1782 940436

www.dwkltd.com
enquiries.uk@dwk.com

PRODUCT SPECIFICATION/CERTIFICATE OF CONFORMITY

Specification	**Product Code**	327150
	Description	30ml Universal Container with Plain label PP Screw Cap
	Material Base	Polypropylene
	Material Cap	Polypropylene
	Colour	Clear
	Sterility	Sterilised by ETO
	Shelf Life	5 Years from Manufacture

Dimension	**Height (mm)**	93.5
	Diameter (mm)	30.5
	Packaging Format	400/Case 50/Bag

Compliance	**Food safe (EU/FDA)**	EU 10/2011 / FDA CFR 21 177
	BSE / TSE free	YES
	REACH Compliant	Compliant according to 1907/2006
	LATEX Free	YES
	CE Marking	Meets the essential requirements of European directive 98/79 EC
	Leak Testing	In Accordance with to BS EN 14254 Annexe D

DWK Life Sciences Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of DWK Life Sciences Limited by:

Signed: K Marshall
Quality Engineer

Date: 23rd November 2020