
CATALOG NO.	362694		362695		362696		362697	
NOM. CAP (mL)	15		15		50		50	
MATERIAL	PP		PP		PP		PP	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.)	0.58	14.7	0.58	14.7	1.07	27.2	1.07	27.2
B (HT W/ CLOSURE)	4.76	120.9	4.76	120.9	4.55	115.6	4.55	115.6
C (HT W/O CLOSURE)	4.71	119.6	4.71	119.6	4.50	114.3	4.50	114.3
D (HEIGHT TO LID)	4.47	113.5	4.47	113.5	4.26	108.2	4.26	108.2
E (DIA.)	0.64	16.3	0.64	16.3	1.10	27.9	1.10	27.9
W.T. (NOM. WALL THK.)	.035	0.89	.035	0.89	.035	0.89	.035	0.89
ANGLE	15°		15°		35°		35°	

CATALOG NO.	PACKAGING	QTY/PKG
362694	BULK PACKED	50
362695	RACK PACKED	50
362696	BULK PACKED	25
362697	RACK PACKED	25



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

DRAWN	MC	APVD PVC	PVC 1/3/05		
CHKD	HRH 1/3/05		1/3/05	PSS 36269X	1204
DOCUMENT					

NUNC **EZ*Flip*** FLIP-TOP TUBES