Cytiva Bio-technology (Hangzhou) Co., LTD
Part of Danaher
1568 Chunjiang East Road,
Tonglu, Zhejiang Province, China
Zip: 311500
T: (86571)6996 5100
F: (86571)6996 5888
Site web: www.Cytiva.com

Key Property Specification

Product ID: Grade GFC

Property	Range	Ideal
GRAMMAGE, g/m^2	50-60	55
THICKNESS, μm@53KPa	225-295	260
POROSITY, s/300mls/1sq in	13-27	20
WET BURST, inches H_2O	≥5.2	8.2