Spoon for foodstuffs, blue





The SteriPlast® spoon for foodstuffs, made of polystyrene, is produced in our class 7 clean room (10,000), individually packaged for disposable use and then sterilised by gamma rays. The spoon for foodstuffs is available in two sizes, 2.5 ml (similar to a teaspoon) and 10 ml (similar to a tablespoon), and is ideal for sampling powders, granulates, pastes and fluids.

- PS, blue
- Individually packaged and sterilised by gamma rays
- Production, assembly and packaging according to clean room class 7
- Complies with EU food requirements and FDA regulations

Spoon for foodstuffs, blue

Length mm	Capacity ml	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.	
127	2.5	100		10		5378-3011	N4
170	10	100		10		5378-3012	N4

Food spoons, long handle, blue



Food spoons, long handle, blue
The long handle makes it easy to reach the bottoms of larger containers as well and allows for sampling of drums, jerrycans and bottles with small opening diameters.

Sampling spoon curved, long handle, blue
The long handle of the plastic spoon makes it possible to empty residue from large containers easily and can also reach into places that are hard to reach such as jerrycans with narrow necks.

- PS, blue
- Individually packaged and sterilised
- Handle length 300 mm
- For single use
- Clean room manufactured (clean room class 7)
- Complies with EU food requirements and FDA regulations

Food spoons, long handle, blue

Total length mm	Capacity ml	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.	
355	5	10		10		5378-5003	N4
375	20	10		10		5378-5007	N4

Sampling spoon curved, long handle, blue

Total length mm	Capacity ml	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.	
360	10	10		10		5378-5011	N4